

Monday, 8 June 1998

1998 IMS Workshop WMG:

**Cryogenics: A New Beginning**


# Setting Up and Calibration of a Cryogenic Test Station

J. M. Pond

Code 6851, Naval Research Laboratory  
Washington, DC 20375


pond@chrisco.nrl.navy.mil





## Outline:

- Introduction
- Overview of a system
- The cryogenic refrigerator
- The vacuum chamber
- Providing a vacuum to insulate the experiment
- Temperature monitoring and control
- Getting an isothermal environment
- Getting microwave signal in and out
- Calibration issues
- Some references and conclusions




## Motivation:

- Provide some useful information to rf and microwave engineers who:
  - have little or no experience with cryogenics
  - have little or no experience with vacuum equipment
  - wish to save themselves some grief by finding out from someone who has been greatly grieved what does and doesn't work
- Initiate and stimulate discussion of techniques for improving cryogenic microwave measurements and testing methods
- Firm belief that, even if HTS hasn't "lived up" to the popular hype, cryogenic electronics is a near term reality:
  - colder is better — higher conductivity, higher mobility, less noise
  - cryocooler and packaging technology is becoming affordable


## Disclaimers:

- Where possible information about vendors and suppliers has been provided. In no case should this be presumed to be an exhaustive list nor a recommendation or endorsement



## System Options:


- Closed-cycle refrigerator
  - Vibration — can be a problem — microphonics
  - Self contained — no transfer of cryogenic fluids — no cryogenic infrastructure required
  - Reliable — based on cryopumps found in vacuum systems
  - Not cheap
  - Ensuring an isothermal environment is more difficult
- Open cycle — need cryogenic fluid (liquid Nitrogen or Helium most common)
  - Nitrogen is cheap and easy to handle but can't get much below 77K
  - Good if you can immerse in LN and are happy with one temperature and can stand condensation on device
  - Liquid Helium requires much more expertise in handling — can get below 4K but is more expensive
  - With either, if you want temperature control over wide range you approach complexity of closed cycle system




## Components of a Closed Cycle System:

- Refrigerator
  - Compressor — where heat is expelled — air or water cooled
  - Expander — the gas expands and cools the “cold finger”
  - Gas lines — often the expander and compressor are separate units with two Helium gas lines connecting them (He plays the role that Freon plays in your household refrigerator)
- Vacuum system — vacuum provides the insulation that allows the DUT to be cooled
  - Pump — various options from roughing pumps to turbo pumps
  - Gauge — allows monitoring of chamber pressure
- Temperature regulation — refrigerator runs unregulated — heaters provide ability to adjust temperature
  - Heaters — resistors
  - Sensors (thermometers) — many options
  - Controller — feedback loop between heater and sensor
- Vacuum chamber — provides thermal isolation of DUT
  - Feed throughs — electrical, mechanical, optical





# Cryogenic Microwave Test Station



Compressor

Vacuum  
Pump


Temperature  
Controller


Vacuum  
Monitor






## The Chamber, Microwave Feedthrough and Temperature Controller






## The Basic Subsystems of a General Purpose Cryogenic Microwave Test System








## Some of the Basic Components




Basic Expander Module



FTCO-12

FTC-K150-18


Basic Vacuum Chamber



Viewport/Optical




Sliding Coax Feedthrough



## Laboratory Closed-cycle Refrigerator:

- Assuming size, weight, and power consumption are of minimal importance in laboratory setting, need to estimate:
  - Minimum operating temperature
  - Maximum heat load at that temperature
    - generated within the DUT/sample
    - due to thermal radiation and thermal conduction
- Many closed-cycle systems for laboratory use have been developed from the cryopump industry:
  - Usually Gifford-McMahon cycle
  - Typical cooling capacities: 10-50 Watts @ 77K and 5-15 Watts @ 20K
- Some vendors:
  - APD Cryogenics ([www.apdcryogenics.com](http://www.apdcryogenics.com))
  - CTI-Cryogenics ([www.ctivacuum.com/enhanced/index.htm](http://www.ctivacuum.com/enhanced/index.htm))
  - CVI ([www.chart-ind.com/cvi/crs/](http://www.chart-ind.com/cvi/crs/))
  - Leybold Cryogenics ([www.leyboldcryogenics.com](http://www.leyboldcryogenics.com))




## Vacuum Chamber:

- Rely on standard commercial parts as much as possible
  - Stainless steel components readily available
    - stainless is the standard for UHV — minimal outgassing
    - stainless is heavy
    - wide range of standard component (building blocks) available to minimize custom work
  - Custom stainless steel vacuum equipment is expensive
    - welding to UHV standards is not trivial
    - brass, aluminum, lucite, etc. are fine for custom/home made parts
- Some vendors:
  - MDC Vacuum Products Corp. ([www.mdc-vacuum.com](http://www.mdc-vacuum.com))
  - Nor-Cal Products ([www.snowcrest.net/norcal/index.html](http://www.snowcrest.net/norcal/index.html))






## Close-up of the Chamber, Cold Stage and Si Diode Sensor









## Top View of Cryogenic Chamber







## Vacuum Chamber — Feedthroughs:

- The ports and feedthroughs are as important as the chamber itself
  - Flexibility
    - have plenty of standard ports using one (maybe 2) flanges
    - brass, aluminum, etc. are fine for custom parts
  - Customizable
    - standard flanges always have blank-off plates available
    - an excellent starting point for custom feedthroughs
  - Expandability — you always need more feedthroughs
    - do not risk integrity of the basic chamber
    - your needs will grow and it is best to be able to keep the basic system geometry intact
- Types of feedthroughs (besides the obvious: coax cable and bias)
  - Mechanical
    - tuning screws on DUT, etc.
  - Optical
    - visual inspection — troubleshooting — show and tell
    - optical response
  - Radome



## Close-up of 6-port Sliding Coax Feedthrough





## Simple Mechanical Feedthrough for Tuning DUT









## A Plexiglass Desiccator Cover for a Radome or for Visual and Mechanical Access






## Vacuum Considerations:

- With decent cooling capacity a rough pump down to ~100 milliTorr is sufficient. The refrigerator will then act as its own cryopump.
- Since the refrigerator is a better vacuum pump than a roughing pump, “backstreaming” of oil into the chamber can be a problem
  - oil can coat the chamber, DUT and cabling — a mess at best
  - solutions:
    - valve off roughing pump — automate by using a pressure gauge/controller ([www.granville.com/275cam.htm](http://www.granville.com/275cam.htm), [www.mdc-vacuum.com](http://www.mdc-vacuum.com))
    - use turbo pump to evacuate chamber
- With a good roughing pump and large refrigerator, virtual leaks, out-gassing and very small leaks are usually not a problem.
- Vent with Nitrogen to eliminate condensation on cold stage/DUT
- Some vacuum pump/system vendors:
  - Alcatel - Balzers - Edwards - Leybold - Sargent Welch - Varian - Precision - Shimadzu - Seiko Seiki - Kinney - Stokes - CTI



## Temperature Control


- Most systems use a refrigerator running at full capacity and use resistive heaters in conjunction with thermometers and a controller to provide an isothermal testing environment.
- Standard practice is to use two thermometers
  - Temperature control – mounted close to the heater
  - Sample temperature – mounted on or as near as possible to the DUT/sample
  - $\Delta T$  of thermometers is a good measure of thermal design
- To provide good control at all temperature:
  - Heater power > cooling power at all temperatures
  - Ohmite “Metal-Mite” aluminum housed axial lead wirewound resistors are a very robust heater (E.J. Cukauskas)
  - “PID” controller provides best performance
  - Autotuning features can be useful
- Cryogenic Temperature Controller manufacturers include
  - Lakeshore Cryotronics ([www.lakeshore.com/](http://www.lakeshore.com/))
  - Scientific Instruments, Inc. ([www.sci-inst.com/](http://www.sci-inst.com/))



## Thermometers/Sensors:


- Types:
  - Diode — fast response, interchangeable, standard calibration
  - Resistor — larger, mounting strain induced errors, many types with both PTC and NTC
  - Capacitor — drift problems and thermal cycling problems, immune to magnetic fields
  - Thermocouple — errors due to temperature gradient of wire
- Characteristics — excellent product and Application Notes are available from vendors — make sure your choice meets your needs:
  - Magnetic field sensitivity
  - Size and Mass
  - Reproducibility and Stability
  - Accuracy/sensitivity/resolution
  - Interchangeability
  - Response time
- Recommendation:  
Unless **H** fields are part of your experiment use 4-wire Si diode






## The Heater Resistors and Control Sensors


(This is the underside and mounts directly on the cold finger)







## Ensuring Good Thermal Contact:

- Minimize number of interfaces that heat must cross.
- Utilize materials such as Oxygen Free High Conductivity (OFHC) copper in as much of the cold finger and DUT mounting platform as possible.
  - not easy to machine
  - particularly important for working much below 77K
  - [www.copperandbrass.com](http://www.copperandbrass.com)
- In all cases it is preferable to directly mount (screws) or clamp the DUT/sample to the cold stage rather than rely on gravity.
- Once the chamber is under vacuum, thermal contact of DUT/sample to cold stage must be ensured. Options include:
  - Thermal grease — Apeizon Type N or equivalent — messy
  - Rubber cement — particularly good for small sample where concern about use of solvents exists
  - Indium foil — more difficult to use — may “stick” to sample and be destroyed in removing — expensive
  - Au foil/indium foil — even more expensive — Au is inert



## Thermal Behavior of Cryogenic Microwave Test Facility Temperature Versus Time and Cooling/Heating Rates






## RF Feedthroughs:

- Balance the issues of attenuation versus heat leak.
- Heat leak not only loads the cooling capacity of the cooler it also is dumping heat directly into the region of the DUT/sample that is being measured leading to higher temperature uncertainty.
- In many laboratory cases, trading off attenuation for better thermal conditions is preferable — the real world is much more challenging.
- Stainless steel coax is popular option compared to Cu


Cable Type	Thermal load (W•cm)	Insertion Loss (dB/cm)
0.085 Cu	14.97	0.022
0.141 SS	1.77	0.028
0.085 SS	0.90	0.048







## Coaxial Cable and Connector Considerations:

- Teflon has high coefficient of thermal expansion
  - gaps can occur at low temperatures
  - if possible use coax with special dielectric with low coefficient of thermal expansion (ISOCORE by Rogers is no longer made)
- Hysteresis in teflon creep can lead to problems:
  - Center conductor contact problems – “captured center pin” connectors may be required especially for long straight cables.
  - Measurement repeatability and, hence, calibration problems
- Any of the standard coax connector series (SMA, K, OS-50) can be used — individual connectors within a series may be suboptimal
- Stainless cables require stainless compatible solders and fluxes. Suggested solutions:
  - Eutector Flux 157 and EutecRod 157 (Eutectic Corporation, 40-40 172nd St., Flushing, New York 11358)
  - Indium Corporation of America ([www.solder.com/index.html](http://www.solder.com/index.html))




## Cold Stage with DUT Mounted Background is 10 Layer MLI






## Measurement and Calibration Issues:

- Everything has some amount of temperature dependence — measurement calibration is an issue
  - If nothing else, the physical length changes
  - IL of coax decreases — electrical length decreases
- Room temperature calibration
  - reduced accuracy
  - easy — same as always
  - some first order correction for IL can be made in some cases
- Cryogenic calibration
  - Standards, particularly loads, have temperature dependence themselves — TRL preferable
  - Tedious — unless a matched network using switches is used, a cryogenic cycle is needed for each standard
  - Instrument drift may be an issue
- Cryogenic calibration
  - Standards, particularly loads, have temperature dependence themselves — TRL preferable




## Measurement Error at Various Temperatures Illustrated by Measurement of an SMA F-F Bulkhead Adapter After Room Temperature Calibration

(using HP 85052D (SOL) 3.5mm economy cal kit)









# Measurement Error at Various Temperatures Illustrated by Measurement of an SMA F-F Bulkhead Adapter After Room Temperature Calibration


(using HP 85052D (SOL) 3.5mm economy cal kit)






## Measurement and Calibration Issues (cont.):

- The previous example was a worst case in many ways:
  - The cable length subject to a thermal gradient was at a maximum (~20 cm)
    - IL can be reduced by ~2X by using 0.141
    - Use shorter cables — 5 cm of 0.141 cable would have a 0.5 Watt thermal load but would reduce IL by ~4X
    - phase error would, similarly, be reduced by ~4X
  - Based on measurements shown in previous graphs, these modifications should result in IL error  $< \sim 0.1$  dB and phase error  $< \sim 0.5$  degrees
- Possible complications:
  - Thermal design of DUT/sample becomes more of an issue since more heat is entering the DUT coaxial connectors.
    - minimize layers and interfaces in DUT package
    - place sample temperature sensor near coax connector
    - more sample temperature sensors to monitor  $\Delta T$  of DUT
  - May need to customize coax feedthroughs to optimize this performance for each DUT



## Measurement and Calibration Issues (cont.):


- Noise measurements can be particularly tricky to the unwary
  - Standard Noise Figure Meters:
    - assume that both ports of the DUT are available
    - are scalar instruments — hence have problems with electrically long low-loss devices
    - cryogenic feedthroughs represent just such a problem — additionally their electrical length and loss are temperature dependent
    - will happily report negative Noise Figures
  - Proceed carefully from first principles
    - even by deembedding from the cascade doesn't really solve the problem due to reflections
    - an isolator can help with the reflection problems and make deembedding more accurate in some cases
  - Listen carefully to the two talks in this Workshop on this important topic
    - low Noise Figure is one of the primary reasons cryogenic operation is so promising
    - it is the measurement area where standard commercial RT techniques are least applicable



## Some Suggested References:

- Cryogenic Engineering, R.B. Scott, D. Van Nostrand Co. Inc., Princeton, New Jersey, 1959
- Experimental Techniques in Low-Temperature Physics, Third Edition, G.B. White, Oxford Science Publications, Clarendon Press, Oxford, 1979
- Cryogenic Process Engineering, K.D. Timmerhaus and T.M. Flynn, Plenum Press, New York, 1989
- Temperature Measurement and Control Catalog, Lakeshore Cryotronics, Inc., Westerville, Ohio
- Web sites:  
    [www.omega.com/techref/itmp.html](http://www.omega.com/techref/itmp.html) — The International Temperature Scale of 1990 (ITS-90)





## Conclusions:

- Performance advantage of cryogenic operation is obvious.
- Cryocoolers for product are increasing in reliability and decreasing in price — particularly if LN temperatures are not needed.
- Cryocoolers for laboratory use are very rugged and reliable.
- Assembling a laboratory system is basically an exercise in assembling off-the-shelf hardware from different fields.
- Custom hardware can be minimized
- Measurement and calibration of a cryogenic system at microwave frequencies, although more difficult than at room temperature, involves basic trade-offs of speed versus accuracy — analogous to issues that affect room temperature measurements as well — just not as familiar to the microwave engineer and there is not, as yet, a solution offered by the major microwave test and measurement industry.